collacone®

Gebrauchsinformation
Instructions for use
Instructions d'utilisation
Instrucciones de uso
Istruzioni per l'uso
Informações de utilização
Használati útmutató
Informacija o korištenju
Kullanım bilgileri
Информация по применению
Πληροφορίες χρήσης

Fabricant:
Manufacturer:

botiss biomaterials GmbH
Hauptstraße 28
15806 Zossen
Germany

Gebrauchsinformation collacone®

DE

Resorbierbares, Lokales Hämostyptikum
Steril

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Sterile

For single use only

Porcine collagen manufactured
from porcine corium

100 g contain a minimum of 96,75 g collagen
of porcine origin

Interaction with other agents

Antiseptics, which release chlorine (e.g. chloramine) as well tannic acid and caustic which change the protein must not be used together with collagen.

Ointments, powder as well silicon preparations can change space at collagenfoil, -sponge and -fleece, and should not be used together with collagen.

Indication areas

collacone® is indicated wherever capillary, venous, small arterial and diffuse seeping bleedings must be stopped, and where conventional means of hemorrhage control are either inadequate or technically difficult and time-consuming.

Indications

collacone® is used as a haemostatic agent following tooth extractions. In patients at risk of developing haemorrhages, the combined use with fibrin sealant is indicated.

- Socket preservation
- Hemostyptic treatment of extraction wounds of patients at risk of bleeding
- Sinus floor elevation
- Periimplant bone defects when placing immediate implants
- Filling the expansion gap in bone splitting
- Filling of bone defects (e.g. after cystectomy)

Dosage and Application

After opening of the package, collacone® is taken out with dry sterile instruments, placed on the cleaned wound and pressed on lightly with swabs. collacone® can be cut to the necessary size with sterile scissors. Immediately after contact with the wet wound surface

collacone® sticks to the wound and forms a gel like bond with blood. With superficial wounds an additional layer of cotton gauze fixed with elastic bandages is recommended.

Contraindications

collacone® is contra-indicated in infected wounds, and must not be used concurrently with methyl methacrylate adhesive, it must not be used for intravascular application.

If you use an absorbable hemostatic agent on or near bony or neural spaces:

- use the minimum amount necessary to achieve hemostasis; and,
- remove as much of the agent as possible after hemostasis is achieved.

This will reduce the likelihood of neural and other soft tissue damage from swelling of the absorbable hemostatic agent, and/or migration and swelling of fragments of the agent.

Side effects

Allergic reactions may occur in isolated cases. Sometimes pain appears after application of a dry preparation to the wound surface. The increase of existing infection is very rare.

Note!

collacone®

- is sterilized by γ -irradiation and must not be resterilized
- must be considered insterile and must be discarded if the package is damaged
- must not be used after the expiry date
- must not be stored above 24°C
- must be protected against moisture
- is a medical device and therefore must be kept from children!

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