dental
bone & tissue
regeneration

botiss
biomaterials

maxresorb®

DE GEBRAUCHSINFORMATION
EN INSTRUCTIONS FOR USE
FR MODE D'EMPLOI
ES PROSPECTO
IT INFORMAZIONI PER L'USO
PT INSTRUÇÕES DE UTILIZAÇÃO
HU HASZNÁLATI ÚTMUTATÓ
HR UPUTE ZA UPORABU
TR PROSPEKTÜS
RU ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ
EL ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ

CE 0483

Fabricant / Manufacturer
botiss biomaterials GmbH
Hauptstraße 28
15806 Zossen
Germany

Rev.: B01MR-02/2015-04

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Please read these Instructions for Use in their entirety before using the product!

The manufacturer is not liable for damages caused by non-observance of these Instructions for Use.

Please note: maxresorb® must only be used by qualified and appropriately trained personnel. Sole responsibility for proper selection of patients, for adequate training, for experience in the use of maxresorb®, for all aspects of the surgical interventions, and for the choice of post-surgical procedures rests entirely with the surgeon.

1. Product description / material
maxresorb® is nominally composed of 60 % hydroxylapatite (HA) and 40 % beta-tricalcium phosphate (β -TCP). These materials have undergone extensive and successful clinical testing over many years. maxresorb® is a safe and biocompatible material that works as osteoconductive scaffolding material to support the ingrowth and fusion of adjoining, vital bone when brought into an osseous environment in accordance with the concept. maxresorb® is resorbed and replaced by natural bone. In many cases, maxresorb® can save patients the trauma of removal of autologous bone. It also offers an alternative to the use of human or animal cadaver bone without the risk of disease transmission.

2. Areas of application

maxresorb® is a synthetic, osteo-conductive bone substitute that is applied in granulate form. maxresorb® is especially intended for the filling of intraosseous, periodontal and root resection defects, for filling of extraction alveolae, for elevation of the maxillary sinus floor and for augmentation or reconstructive treatment of the alveolar ridge.

maxresorb® can be used in conjunction with products intended for Guided Bone Regeneration (GBR). The granulate can be applied directly into the defect, or it can be mixed with sterile saline solution, blood, bone marrow aspirate or autogenous bone before implantation. maxresorb® works as temporary scaffolding and is not intended to provide structural support during the healing process. The bone replacement material is biocompatible, radiopaque, and is resorbed in a controlled manner over the course of time.

3. Contraindications

maxresorb® is only intended for the specified uses and excludes any other possibilities for use. maxresorb® must not be used in the presence of a contraindication. maxresorb® is contraindicated as a structural support and must not be used for fixation of titanium anchors.

Additional contraindications:

- Existing acute or chronic infections, particularly at the implantation site
- Severe vascular or neurological diseases
- Non-controlled diabetes
- Severe degenerative diseases
- Hypercalcaemia, abnormal calcium metabolism
- Inflammatory bone diseases
- Malignant tumours
- Severely impaired renal function

4. Precautions

In order to achieve appropriate bone regeneration it is important to ensure

that there is maximum contact between the existing bone and the bone replacement material. maxresorb® must not be used for direct bone anchoring. Titanium anchors must not be implanted in the augmentation area before approx. 6 months have elapsed after application of maxresorb®.

5. Possible complications

A successful result cannot be achieved in every case. A secondary operation for removal or replacement of the bone replacement material utilised may become necessary due to a surgical error, certain medical circumstances, or the failure of the bone replacement material. Possible adverse effects may include, but not be limited to, the following:

- Fracture or compression of the bone replacement material, with or without formation of debris due to a strong application of force;
- Bone deformation/non-viable

- bone at the implantation site;
- Post-surgical and other complications that may occur during a surgical intervention;
- Allergic reaction to the product, although no adverse reactions have been reported to date.

6. Storage and durability instructions

Storage should be in a dry location at room temperature. The expiry date is printed on the label. Do not use maxresorb® after the expiry date. maxresorb® is environmentally friendly.

Special disposal is not necessary.

7. Application

Step 1: Uncover the defect with a mucoperiosteal flap and carefully remove the existing granulation tissue. Open the outer and inner pouch of the maxresorb® packaging. Open the sterile primary packaging.

Step 2: The bone replacement material may be first mixed with sterile saline solution, blood, bone marrow aspirate or autogenous bone in accordance with customary surgical procedures and then applied to the corresponding area. An excessive use of force can result in compression of the trabecular structure. Overfilling of the defect must be avoided. The bone replacement material must only be brought into direct contact with well-vascularised bone. Cortical bones should be mechanically perforated.

Step 3: After implantation, the surgical field should be secured to prevent migration of the implant. The mucoperiosteal flap created should be allowed to heal by first intention through placement of a single suture. The wound can be covered with a bandage for a period of one to two weeks. Approx. 6 months will be

required before the bone satisfactorily colonise the bone replacement material and the titanium anchors can be inserted.

8. Warnings

The contents are in double-sealed packaging and are sterile. If the packaging is open or damaged, the bone replacement material must not be used. The product must not be used after its expiry date. It is recommended that the inner bag should not be opened until a few minutes before implantation. maxresorb® is only intended for single use. It must not be resterilised or reused. The granulate must be safely affixed to prevent migration. This must be taken into account when selecting an augmentation procedure. maxresorb® is radiopaque and could possibly hide areas above or below the bone replacement material on an X-ray.

9. Further instructions

These Instructions for Use satisfy the current state of technology and our own experience. The product must only be used for the indications described under section 2. The user is personally responsible for use of the product. The manufacturer is not liable for defective results since it is not involved in the surgical workmanship.

Any claims for damages that may nevertheless arise shall relate exclusively to the merchandise value of our products.

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