dental
bone & tissue
regeneration

botiss
biomaterials

maxresorb[®] inject

DE GEBRAUCHSINFORMATION
EN INSTRUCTIONS FOR USE
FR MODE D'EMPLOI
ES PROSPECTO
IT INFORMAZIONI PER L'USO
PT INSTRUÇÕES DE UTILIZAÇÃO
HU HASZNÁLATI ÚTMUTATÓ
HR UPUTE ZA UPORABU
TR PROSPEKTÜS
RU ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ
EL ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ

CE 0483

Fabricant / Manufacturer
botiss biomaterials GmbH
Hauptstraße 28
15806 Zossen
Germany

Rev.: B01MI-02/2015-04

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Description

maxresorb® inject is a safe and fully biocompatible material that is used as an osteoconductive agent to support the ingrowth and fusion of adjoining, viable bone material, when used in an osseous environment. Due to a structure that is fundamentally identical to that of the mineral phase of human bones, maxresorb® inject is fully resorbed by natural bone material and replaced by it.

The hydroxylapatite (HA) component in the paste is very quickly resorbed, while the hydroxylapatite / tricalcium phosphate granule is resorbed more slowly and thus supports the bone with a mineral matrix for a longer time. maxresorb® inject can often save patients the trauma of an autologous bone transplantation and is a sterile, fully synthetic alternative to human or animal bone material and, in contrast, cannot transmit any diseases.

Composition

maxresorb® inject is a combination of 83.5 % HA paste and 16.5 % bone transplantation granules (0.5–1.0 mm). The HA paste consists of 38 % hydroxylapatite in water, the porous granules consisting of 60 % hydroxylapatite and 40 % tricalcium phosphate. Hydroxylapatite and tricalcium phosphate have been extensively and successfully tested in clinical studies over many years.

Areas of application

maxresorb® inject is a synthetic, osteoconductive bone paste for the filling and supporting of adjacent, viable bone material in cases of bone defects that are not critical to the stability of the bone structure. These defects can be found in the extremities, the pelvis, the dental and peridental areas, or the jaw. maxresorb® inject can be used to fill bone defects as a result of fractures, after the removal of

benign tumours or cysts, or to fill transplantation areas. maxresorb[®] inject can also be used for defects of the alveolar ridge and the alveolar wall, including extraction alveolae, and to fill bone defects. maxresorb[®] inject is carefully filled into the defect by hand and, if necessary, moulded and fitted more accurately into the defect. If needed, it can be mixed with bone marrow aspirate before implantation. maxresorb[®] inject acts as a temporary support and is conserved as a soft paste. It does not offer any load support whatsoever during the healing process. The implant is biocompatible, radio-opaque, and is resorbed in a controlled manner.

Contraindications

maxresorb[®] inject must only be used for the specified areas of application. Do not use maxresorb[®] inject if one of the contraindications is present. maxresorb[®] inject must

not be used as a load support for the skeletal system and **MUST NOT BE USED IF THE IMPLANTATION SITE IS UNSTABLE OR NOT FIRMLY AFFIXED, OR IS BEING USED AS A SCREW FIXATION SITE.**

Other contraindications include:

- Acute or chronic infections, especially in the implantation site
- Severe vascular or neurological diseases
- Inadequately vascularised implantation site
- Untreated diabetes
- Severe degenerative diseases
- Diseases of the connective tissue
- Hypercalcaemia, abnormal calcium metabolism
- Inflammatory bone diseases
- Malignant tumours
- Markedly impaired renal function
- Open epiphyseal plates in paediatric patients

Precautions

maxresorb® inject is only intended for professional use by surgeons who have knowledge of bone restoration and bone replacement techniques and who are qualified to perform them. maxresorb® inject is not intended for use in bones providing load support. The implantation site must be mechanically safeguarded using a rigid fixation device or a cage in order to provide a static and load-free environment for the implantation. maxresorb® inject must not be used to affix bone screws. Contact between existing bone and the implant must be maximised in order to ensure correct bone regeneration. maxresorb® inject must only be used with fresh, vital bone material so that rapid revascularisation can be guaranteed. It is important to work with the highest level of care in order to prevent embolisms. For this reason, high-pressure application of

maxresorb® inject in tightly confined defects with venous or arterial access must be avoided. The effect of maxresorb® inject on patients with the following diseases is unknown.

- Chronic infections
- Metabolic bone diseases
- Radiation bone therapy
- Cardiovascular diseases

The effects of maxresorb® inject on paediatric patients and the consequences of mixing it with other substances is unknown. At the time of publication, however, no special precautionary measures have been identified.

Possible complications

A successful result cannot always be guaranteed. A second operation to remove or replace an implant could become necessary due to a surgical error, a particular disease, or equipment failure.

Adverse events

Possible adverse events include, among others:

- Displacement of maxresorb[®] inject due to a load on the bone
- Bone deformation/non-viable bone at the implantation site
- Post-surgical complications such as infections, haematomas, oedemas, swelling, accumulation of fluids, thinning of tissue, and other complications caused by surgery
- Allergic reaction to the product

Warnings

The product is gamma-sterilised and packaged in sealed double-bags. **DO NOT USE IF THE BAG IS OPEN OR PUNCTURED, IF THE PRODUCT AND/OR THE APPLICATOR APPEAR TO BE DAMAGED, OR THE PASTE HAS DRIED OUT.** Check the date of expiry before use and **DO NOT USE**

THE PRODUCT AFTER the expiry date. Do not open the inner bag until a few minutes before implantation. maxresorb[®] inject is only intended for **SINGLE USE**. Do not attempt to resterilise or reuse the product. This can result in cross-infections. maxresorb[®] inject should only be applied using techniques with which the implant can be adequately contained. maxresorb[®] inject is opaque to X-rays and could possibly hide areas above or below the implant on an X-ray.

Application

Step 1: Open both the outer and the inner bags. Remove the applicator's seal cap. maxresorb[®] inject is supplied as a ready-to-use paste in the applicator.

Step 2: Implant the amount of maxresorb[®] inject that is required for the defect. The maximum dosage of maxresorb[®] inject is 20 g per defect. The defect should be filled completely and the paste should be moulded so as to adapt to the natural contours of the bone. The product should make direct contact with all the surfaces of the defect. If the bone to be treated is under load, it should be used together with a rigid fixation device.

Step 3: Protect the implantation site from movement or displacement of the implant. In the event of excessive accumulation of fluids, the area surrounding the implant may be cauterised or suctioned in order to thus reduce haemorrhaging.

Storage

Store maxresorb[®] inject between 5°C and 30°C. Avoid direct sunshine or contact with heat sources.

Shelf life and disposal

The expiry date can be found on the label. Do not use maxresorb[®] inject after the expiry date. maxresorb[®] inject is environmentally friendly and does not require a special disposal method. The used applicator should be disposed of as clinical waste.

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